

Alcohol and Drug Testing

HK Drugcube 10+








The HK DrugCube 10+ is designed to give a preliminary onsite result in just 5-10 minutes. It is used to detect active drugs present in oral fluid and therefore indicates recent use. It also has an added adequacy indicator to show when sufficient sample has been collected.








Non-negative results can be sent to our UKAS accredited laboratory for confirmation analysis using our Chain-of-Custody collection kits (please ask for details).

Hampton Knight can provide competency training, testing protocols and consent forms to ensure your testing process is robust and legally defensible.

The **HK DrugCube 10+** provides an on-site qualitative test for the presumptive presence of the following drug types:

 Alcohol	 Ketamine
 Amphetamine	 Methadone
 Barbituates	 Methamphetamine
 Benzodiazepines	 Opiates
 Cannabis (THC)	 Propoxyphene
 Cocaine	

Key Features

-  10 Minute test result
-  No need to dip reagent strips or cards
-  No need to drip specimen with pipette
-  Each drug type test strip has built in procedural control
-  Built in adequacy indicator to show when sufficient sample has been collected

For further information or to order this product please call 01827 65999 quoting product code: ORT011