

HK Exxon 9










The **HK Exxon 9** is a **Split key drug screening cup** which has a **valve built in with a key to activate** when you are ready to run the test. It is a **combined collection cup and testing device** that **delivers instant results with unrivalled ease and convenience, built to ExxonMobil 2015 specifications.**







Non-negative results can be sent to our UKAS accredited laboratory for confirmation analysis using our Chain-of-Custody collection kits (please ask for details).

Hampton Knight can provide competency training, testing protocols and consent forms to ensure your testing process is robust and legally defensible.

The **HK Exxon 9** panel provides an on-site qualitative test for the presumptive presence of the following drug types:

 Amphetamine (500 ng/ml)	 Ecstasy (MDMA) (500 ng/ml)
 Barbiturates (300 ng/ml)	 Methadone (300 ng/ml)
 Benzodiazepines (300 ng/ml)	 Opiates (300 ng/ml)
 Cannabis (THC) (20 ng/ml)	 Oxycodone (100 ng/ml)
 Cocaine (150 ng/ml)	

Key Features

-  5 Minute test result
-  Uses ExxonMobil 2015 cut off levels
-  No need to dip reagent strips or cards
-  No need to drip specimen with pipette

For further information or to order this product please call 01827 65999 quoting product code: HKDC2016